


#### Quality CPR Feedback:


- Compression depth
- Compression rate
- Compression fraction
- Chest complete release
- Correct hand placement
- Jaw thrust and head tilt control
- Appropriate ventilation volume




## Quality CPR Training

# CPRLilly PRO+

The new **Quality CPR manikin** enables instructors and healthcare providers to measure, monitor and analyze the CPR performance of up to **10 trainees at a time**. CPRLilly PRO+ helps instructors increase effectiveness and efficiency in their CPR courses by connecting with the new **CPRLilly App** on tablet to follow the CPR performance and deliver objective feedback leading to Quality CPR training.



Instructor will at any time be able to access a detailed view on each trainees CPR-performance by a simple touch



Trainees will follow the CPR parameters on a simplified feedback view and improve their performance in real time

### Key Training Features:

- Quality CPR feedback with debriefing history
- Mouth-to-mouth & mouth-to-nose ventilation
- Ventilation with face shield, pocket mask and BVM
- Chest rise with proper ventilation
- Control of proper airway desobstruction
- Carotid pulses
- AED pads placement training capable
- Free App available on google play and iOS
- Constant analysis of CPR parameters for up to 10 manikins simultaneously
- Individual review possibility in group training for more efficient training
- Adaptable to different course concepts
- Simplified student feedback view available with the CPRLilly App

### Solid manikin for effective training:

- Made in Germany of extremely durable, latex free material
- Face masks and lungs bags can be changed easily without tools
- Easy to clean with hygienic airway system
- Runs quietly with efficient power usage and rechargeable battery
- Stable WiFi connection and intuitive App interface
- Practical carrying bag with training mat for comfortable practise
- Dark skin manikin available

### Delivery content:

1 CPRLilly PRO+ training manikin, 1 Tablet (depends on the model), 1 charger and battery, 4 lower face masks (incl. teeth), 25 throat bags, 12 AED-Pad adapters, 1 carrying bag with training mat

Upgrade now and equip your courses with direct feedback CPR training manikins to comply with AHA's guidelines. CPRLilly PRO+ meets the requirements for direct feedback and Quality CPR of the AHA guidelines:

“Instrumented directive feedback device measures compression rate, depth, hand position, recoil and chest compression fraction and provides real-time audio or visual feedback (or both) on these critical CPR skills.”

**Manikin with Tablet 1022071**

**Manikin without Tablet 102237**



9001134\_9/2019



MORE INFORMATION AVAILABLE ONLINE AT [3BScientific.com](http://3BScientific.com)



**3B Scientific GmbH** • Ludwig-Erhard-Straße 20 • 20459 Hamburg • Germany • [3bscientific.com](http://3bscientific.com)

**Phone:** +49 40 73966-0 • **Fax:** +49 40 73966-100 • **E-mail:** [info@3bscientific.com](mailto:info@3bscientific.com)